




Lightweight and Durable Frame




Comfortable Knee Padding



Height Adjustable Knee Pad



Height Adjustable Handle



8" Wheels

The Premium Knee Walker is your key to comfort and mobility when recovering from a below knee injury or surgery. This product is designed to suit your needs with its lightweight, durable frame, comfortable knee padding, large wheels for easy travel over uneven surfaces, and a bag to carry your belongings with you as you travel.

Premium Knee Walker

KA337

Features

- Lightweight and durable frame, easy to manoeuvre through the recovery period
- Contoured knee pad, designed for maximum comfort
- Large wheels, perfect for traveling smoothly over uneven surfaces
- Adjustable handle height and knee pad height to suit the user's needs
- Convenient storage bag
- Non-weight bearing support
- A great addition to your hire equipment pool

Specifications

Overall length	800mm
Overall width	420mm
Overall height	790mm - 1030mm
Handle width	405mm
Knee pad height	430mm - 590mm
Knee pad dimensions	350mm x 150mm x 80mm
Wheels	8 inches
Product weight	10.9kg
Maxi bag load	5kg
Max user weight	136kg

Please contact customerservice@kcare.com.au for more information.
www.kcare.com.au | 1300 783 783