

Afiscooter

*The top of the line mid-size
Multi-purpose scooter*



The **Afiscooter C** with its unique design made to maximize the protection and safety of the driver and its excellent multifunctionality **is the perfect mobility scooter for those who wish to maintain an active, dynamic and versatile lifestyle.**

The 3 wheel model is supremely maneuverable and the 4 wheel model is ultra-solid, offering smooth and comfortable ride along with improved stability.

You will be travelling in outstanding style and comfort, Wherever you wish.

Highlights

- A robust design to maximize the protection and safety of the driver.
- A delta tiller with a LCD display, left and right side finger touch controls to accommodate any driver.
- Full front and rear suspension system for ultimate comfort while driving (for both indoor and outdoor use).
- A luxurious, orthopedic and remarkably conformable seat, provides ease boarding and comfort while driving.
- Excellent storage space with a lockable front compartment and an option for an additional lockable rear box.
- A powerful, long-lasting, energy efficient led headlight.
- An option for a hard-top canopy to ensure comfortable ride in all weather conditions.



Pleasurable
driving experience



Design for
safety



Ease of
operation



Top quality
& Reliability



Afikim electric vehicles is a leader in developing and manufacturing premium, electrically-powered mobility scooters for more than 40 years. The well-known Afiscooter product line is designed with driver safety in mind first and foremost, alongside attention to ergonomics, comfort and reliability of the vehicle. Afikim electric vehicles has established a worldwide distribution network that offers local technical support and follows the highest standards of service.

Our company meets highest standards of quality and European ISO 9001, CE and EN 12184 certifications, as well as US safety standards with FDA approval.

Your independence
is not a matter
of weather



Features



Specification

		Afiscooter C3	Afiscooter C4
Overall Length		130cm (51")	137cm (54")
Overall Width	Standard Wheels	67cm (26")	
Overall Height	Folded	80cm (31")	
	(-) Canopy	128cm (50")	
	(+) Canopy	163cm (65")	
Ground Clearance		10cm (4")	
Weight	(-) Batteries	75kg (165lb)	95kg (209lb)
	(+) Canopy	+12Kg (26lb)	
Weight Capacity	Incl. Driver	150kg (330lb)	180kg (400lb)
Travel Range	(+) Batteries	Up to 50km (31m)	
Maximum Speed		10-15Kph (6.2-9.3mph)*	
Max Climbing Angle		9° (16%)	
Turning Radius		128cm (50")	170cm (67")
Seating Width		46cm (18")	
Tire Type		Pneumatic (1.7 bar)	
Tire Size	Front Tyres	4.10/3.50-6"	
	Rear Tyres	4.10/3.50-6"	
Battery Type	two	12V 85Ah, Max.	
Power Unit	Rear wheel drive	24V DC motor 950W	
Controller Type		SD 120A	
Charger Type		8 Amp	
Options/Accessories		Canopy Rear lock box Wide seats 20" Crutch Holder	

* Speed limit varies between the different countries according to the local regulation of each country.

** U.S. model | Tested by manufacturer, based on actual parameters | The above specifications subject to change without notice