




Lightweight and high strength frame



Bag under seat for storage

# Bass Magnesium Walker

## Features

- Lightweight and high strength Magnesium Alloy frame
- Height adjustable handle with indication numbers
- Height adjustable seat with multiple height options and three factory presets
- Large 8" wheels
- Depth adjustable backrest
- Ergonomically designed handle grips
- Easily side foldable frame
- Bag under seat for storage
- Removable seat pad for easy cleaning
- Tested to ISO11199-2:2021
- 135kg MUW
- Cup holder, cane clip and cane bucket optional

## Specifications

	<b>KA380MGS</b>	<b>KA380MGM</b>	<b>KA380MGL</b>
Product Code	<b>S</b>	<b>M</b>	<b>L</b>
Seat Height	537mm	560mm	605mm
Product Net Weight without Bag & Backrest	6.8kg		
Handle Height	850-990mm		
Overall Width	625mm		
Overall Length	685mm		
Minimum Turning Radius	750mm		
Handle Width (centre to centre)	480mm		
Seat Size	W460 x D230mm		
Folded Dimension	L685 x W250 x H850mm		
Wheels	8"		
Bag Load Capacity	5kg		
Max User Weight	135kg		
Main Frame Material	Magnesium Alloy		

Please contact [customerservice@kcare.com.au](mailto:customerservice@kcare.com.au) for more information.  
[www.kcare.com.au](http://www.kcare.com.au) | 1300 783 783