



Lightweight and high strength frame



Bag under seat for storage

Flinders Carbon Fibre Walker

Features

- Lightweight and high strength Carbon Fibre frame
- Height adjustable handle with indication numbers
- Height adjustable seat with multiple height options and three factory presets
- Large 8" wheels
- Depth adjustable backrest
- Ergonomically designed handle grips
- Easily side foldable frame
- Bag under seat for storage
- Removable seat pad for easy cleaning
- Tested to ISO11199-2:2021
- 135kg MUW
- Cup holder, cane clip and cane bucket optional

Specifications

Product Code	KA380CFS	KA380CFM	KA380CFL
Seat Height	S 515mm	M 560mm	L 605mm
Product Net Weight without Bag & Backrest	5.7kg		
Handle Height	850-990mm		
Overall Width	600mm		
Overall Length	685mm		
Minimum Turning Radius	750mm		
Handle Width (centre to centre)	480mm		
Seat Size	W460 x D230mm		
Folded Dimension	L685 x W250 x H850mm		
Wheels	8"		
Bag Load Capacity	5kg		
Max User Weight	135kg		
Main Frame Material	Carbon Fibre		

Please contact customerservice@kcare.com.au for more information.
www.kcare.com.au | 1300 783 783