



Extra heavy-duty construction



Removable splash guard



Width adjustment pins



Over Toilet Frame Heavy Duty

Width Adjustable 3 in 1

KA450HDWA

This heavy duty over toilet frame has been designed around comfort and stability, for your peace of mind. The height adjustable legs and width adjustable frame ensure that the product fits anywhere and suits your unique needs. The armrests make it both easy and comfortable for transfer, and the anti-slip feet ensure your stability while seated or standing. The seat material and removable splash guard make cleaning simple, and effortless. This versatile product can also be used as a shower chair or bedside commode if the bucket is used.

Features

- Ultra-light and extra heavy-duty construction
- Height adjustable legs
- Width adjustable frame to meet the needs of the occupant
- Comfortable padded seat and armrests
- Removable splash guard for easy cleaning
- Easily fits over most toilets
- Anti-slip tips for stability
- Easy tool-free assembly
- Can also be used as shower chair or bedside commode
- Bucket with lid included for use as a bedside commode
- Toilet roll holder optional
- Maximum User Weight 300kg
- Product tested to AS/NZS 3973:2009 (R2016) Shower/ Toilet Chairs (mobile and static for product safety)

Specifications

Product Net Weight	4.6kg
Overall Height	677-777mm
Overall Width	600-750mm
Overall Depth	465-480mm
Seat Height	500-600mm
Seat Size	460mmW x 415mmD
Width between Side Frame	500-650mm
Maximum User Weight	300kg

Please contact customerservice@kcare.com.au for more information.
www.kcare.com.au | 1300 783 783