




Yoga

Yoga

Features

- Strong light-weight alloy aluminum frame designed for easy folding and carrying
- Small foldable scooter for easy storage and transportation
- Lightweight 24 volt Li-ion battery with quick charge 2A off board charger
- Height and angle adjustable steering tiller to fit individual users
- Ideal for both indoor and outdoor use
- Patented Twin-front-wheels w/spring suspension system offers excellent stability
- Armrest as standard.
- Comfortable padded seat and backrest
- Optional battery charging station, for charging the battery away from the scooter.
- Built in speed reduction switch which is activated when turning offering exceptional safety.

Dimensions cm/in



Overall Dimension – Unfolded : A x B x C
98 x 43.5 x 102.5cm/
38.6 x 17 x 40.3"

Overall Dimension – Folded : A x B x C
52 x 43.5 x 64cm/
20.5 x 17 x 25.2"

Seat Dimension : D x E x F
40.5 x 35 x 31.5cm/
16 x 13.8 x 12.4"

Seat Height from Ground : G 59.5cm/23.4"

Seat Height from Deck : H 39cm/15.4"

Wheel Base : I 78.5cm/31"

Weight kg/lb

battery/ea	total with battery, w/o armrest	capacity
1.7/	25.7/56.7	113/250

Australian Standards Compliant AS/NZS 3695.2:2013

merits
AUSTRALIA

Merits Australia is committed to constantly improve and develop products, we reserve the right to alter specifications without prior notice. Pictures for illustration purposes only.



Side



Side



Hard Case Open



Hard Case



Height adjustable tiller



Folded

Technical Specs

Max Speed	6kph/4mph
Range	18km/11.25mi
Turn Radius	96cm/37.8"
Ground clearance	9.5cm/3.7"
Motor	DC 24V, 80W (Rating)/ 340W (Max)
Controller	PG S-Drive 45A
Battery	Li-Ion Battery 24V/11.5A/H
Charger	2A off-board
Gradient	6°
Front Wheel	6" (153x40) foam filled tire
Rear Wheel	7" (178x45) foam filled tire
Brake	Electromagnetic brake

KCare
Healthcare Solutions