



Lightweight aluminium design

The Economy Bath Transfer Bench is designed to make your bath transfers safer and more convenient. Enjoy the comfort of a supportive seat and a removable backrest. Customise the height to your preference with height adjustable legs that offer smaller increments of adjustment for precise control. Suction cup feet provide security in the bathtub, while anti-slip tips ensure improved stability.

## Economy Bath Transfer Bench

KA224

### Features

- Lightweight aluminium design
- Adjustable seat height
- Removable backrest can be used on either side of the seat
- Moulded plastic seat for easy cleaning
- Suction cup feet for security in the bathtub
- Anti-slip tips for stability
- Easy tool-free assembly

### Specifications

#### KA224

Product Net Weight	4.8kg
Overall Height	800 - 900mm
Overall Width	760mm
Overall Depth	540mm - 570mm
Armrest Height	650 - 750mm
Seat Height	490 - 590mm
Seat Width	670mm
Seat Depth	410mm
Maximum User Weight	136kg

Please contact [customerservice@kcare.com.au](mailto:customerservice@kcare.com.au) for more information.  
[www.kcare.com.au](http://www.kcare.com.au) | 1300 783 783